

Culturing notes from meeting with Paul 3-27-2014

Sydney E. Everhart & Brian J. Knaus

February 27, 2018

Contents

1	Background	2
2	Media	2
2.1	Condensation Management	2
3	Cleaning Cultures	3
4	Short-term Storage	4
5	Long-term Storage	4

List of Tables

1	recipe for cornmeal extract ¹	2
2	recipe for CMA ²	2

1 Background

These are the protocols that we learned from Paul Reeser.

2 Media

Paul's preference is to use cornmeal agar (CMA) for all species of *Phytophthora*. Over the years, however, he's found that commercial formulations of CMA (both Difco and BBL) produce media that is both variable in content and often has negative effects on growth. They now make their own CMA using cornmeal from the Co-Op (Tables 1 and 2). As a part of Paul's protocol he has included using a refractometer to standardize specific gravity of the extract. It was in this process that Paul found cornmeal purchased from places other than the Co-Op (e.g., Winco), yielded lower solid content (chunkies).

Table 1 presents a recipe for cornmeal extract.

Table 2 presents a recipe for CMA.

Table 1: recipe for cornmeal extract¹

?g	Cornmeal
??? mL	water

¹Sydney does not understand what a caption is.

Table 2: recipe for CMA²

?mL	cornmeal extract
?g	agar
?mL	water
??min	sterilization

¹Brian does not understand what a title is.

2.1 Condensation Management

Paul indicated that condensation within Petri plates is a very important step in maintaining *Phytophthora* cultures and should be prevented for a few reasons. Condensation can be conducive to contamination and makes it difficult to identify bacterial contamination. There are several tips that Paul gave us to minimize condensation:

1. Wait until the media has cooled sufficiently before pouring into Petri plates
2. Stack plates high in the hood to minimize condensation
3. Store plates at room temperature in crisper boxes to fully eliminate condensation

- 24 4. Store plates in plastic Petri sleeves at room temperature (never store in
25 refrigerator)
- 26 5. If there's any remaining condensation in sleeve or on plates, allow them
27 to dry in the hood before plating
- 28 6. Petri plates should not be used until they are at least one week old and
29 old media is not a concern (up to a year old can still be used, even with
30 antibiotics amended as long as they are kept in the dark)

31 The benefit of preventing condensation within Petri plates is that it allows
32 para-filming of all plates, which is required in the Hansen lab for all Petri plates.
33 Para-filming helps slow the spread of contaminants that may arise and also helps
34 minimize dehydration of the Petri plates.

35 3 Cleaning Cultures

36 Cleaning cultures is a step that Paul performs for all new isolates sent to him
37 and for any new isolates that are from fresh material. Cleaning cultures is the
38 following steps, which can be done outside the hood unless noted otherwise:

- 39 1. Transfer isolate to dilute CARP (cornmeal, ampicillin, rifampicin and pi-
40 maricin)
- 41 2. Visual inspection of original plug using a dissecting microscope to see if
42 there is any sign of bacteria or other fungal contaminant
- 43 3. Hyphal tipping, performed using a dissecting microscope and needle
- 44 4. Transfer to new CARP plate
- 45 5. After growth, perform Visual inspection for contaminants using a dissect-
46 ing microscope and if any are found, return to step 3
- 47 6. If culture found to be clean, plate on CMA in hood

48 Performing hyphal tipping and inspecting for any signs of contamination
49 was a step that Paul felt was very important for generating clean cultures. His
50 reasoning was that when grown on antibiotic and fungicide amended media,
51 the *Phytophthora* should be able to grow more rapidly than a contaminant. So
52 taking this tissue from the extreme edge of the growing hypha should ensure
53 a transfer of pure culture. Similarly, Paul noted that if a culture was indeed
54 contaminated that most likely places to see signs of contamination is at the
55 site of the original transfer. If a bacterial contamination is suspected and you
56 want to determine if the substance is indeed bacterial, Paul suggests streaking
57 a sample of the material onto nutrient agar.

58 Note that Botran may be used as a non-carcinogenic substitute fungicide
59 for PCNB in PARP.

60 4 Short-term Storage

61 Short term storage of cultures involves the regular transfer of cultures to condensation-
62 free CMA media on a semi-weekly basis. Cultures from the previous week can be
63 harvested by excising a portion of the colony near the growing edge of the colony.
64 This agar plug, containing both mycelia and media can then be transferred to a
65 new plate containing fresh media. Sterile technique should be practiced during
66 this process.

67 5 Long-term Storage

68 Long term storage is accomplished in duplicate, but non-identical, vials. One
69 vial contains deionized water, the second contains 2-3 ‘chips’ of hemp seed (can
70 be purchased at the Co-op). These vials need to be autoclaved for around
71 90 minutes prior to use. Paul uses a special vial which has an indicator in it
72 to ensure that sufficient heat is achieved during sterilization. Also, the hemp
73 seed containing vials should be used within the first three months. If not used
74 before three months they appear to develop something which is detrimental to
75 growth. Three to four agar plugs containing mycelium can be placed in the
76 water containing vial. One agar plug can be stored in the hemp seed containing
77 vial. These vials can be stored at room temperature, in the dark. Cultures can
78 be stored like this ideally for about 2-3 years, but may remain viable for up to
79 15 years.

80 Ken Johnson uses a similar strategy for *P. infestans*, but instead of using
81 hemp seed he uses a few rye seeds.